



TechniCare
Improving Healthcare



Salubrious
Good Medical



Bravo compact bed

A perfect storable, transportable and full functioning hospital bed.

Patented, application number - 2450416-9



Introduction

To improve healthcare around the world Salubrious and TechniCare innovate, develop and manufacture medical technology products. Based on input from nurses, doctors and other healthcare professionals we have developed a fully functional “field hospital bed”, with a strong focus on optimal storage and transportation.

Salubrious is responsible for sales and service of the product and TechniCare is responsible for development and manufacturing.



“We were searching for the optimal compact bed — one that offers efficient storage, is easy to transport, and can be set up without any tools.”





Optimal storage

Storage is one of the most important aspects of the bed.

Since the beds are often not in use, they need to be easy to store.

They can be folded into compact units, allowing up to 7 beds to fit on just 2 pallets.

In total, you can store up to 70 beds in a single 40-foot container.



Easy to install

On-site, it's essential that the beds can be installed easily without the use of any tools.

We have ensured a simple and fast installation process.

The bed can be assembled in approximately 30 seconds—no tools, screwdrivers, or loose components required.

All you need to do is place the upper frame onto the underframe, lock it in place, unfold the upper frame, and you're ready to go.

Technical details

We want to ensure that the beds are robust and still fully functional even after 10 years of use. That's why we choose only high-quality, durable materials. The bed is constructed from stainless steel, features brass bearings, and includes premium components throughout. The wheels are from TENTE and are suitable for outdoor use. In addition, we ensure that maintenance can be carried out on-site using dedicated service packages. We provide a complete list of all components, and every part is easy to replace by any qualified technician.



Specifications

Mattress Surface	2000x800mm
Height adjustment	450mm-900mm
Total weight of the bed	85 KG
Total weight underframe	appr. 40 KG
Total weight upper frame	appr. 45 KG
Total Patient Load	180 KG
Safety working Load	200 KG
Hydraulic actuator	MK5 Power Cable
Bearing Materials	Bronze and POM-C Plastics, chemical resistant.
Frame Material	Steel powder coated, aluminum anodized. Wear parts made from stainless steel 304.
Cover material	Aluminum anodized
Stackable	7 Pcs high on trolley
Adjustment backrest	Lockable gas spring, stepless.
Adjustment CPR	0-85 degrees, lockable gas spring, manual operation by handle bar.
Wheels	TENTE Levina care wheels 125MM
Regulations	MDR 2017/745 Compliant Class ISO13485:2016 ISO14971:2019

Accessories

In addition to our standard configuration, we offer a range of accessories to help you create the ideal bed for your specific needs.

These include various options, from DIN rails and IV poles to the Trendelenburg position. With these accessories, you can configure the bed for both ICU and daycare settings.

Accessories

Headboard	On head and foot section
IV Poles	2 Pieces
Din-rail	On head and foot section
Central break system	Including central brake system and driving position
Adjustment leg section	Fowler position, 6
(Anti) Trendelenburg	Adjustable mechanically
Matress	Anti slip





Improving Healthcare

TechniCare and Salubrious aim to improve healthcare by providing technological innovation.

We deliver high-quality products developed in collaboration with end users. Salubrious is a Swedish manufacturer and distributor of medical devices.

TechniCare is a Dutch company with their own development and production facilities. For more information:

www.technicare.nl & www.salubrious.se



Nikkelstraat 24
2984 AM Ridderkerk

T +31 85 0870 421
info@technicare.nl

www.technicare.nl



Skiffervägen 2G
443 61 Stenkullen

T +46 31 730 58 00
info@salubrious.se

www.salubrious.se